


# RetroView Colonoscopes EC34-i10T

Exceptional control meets enhanced image quality



# To facilitate detection and resection of even hard to reach lesions



### **Optimal therapeutic approach**

Eased therapeutic access to colonic lesions for improved resection procedures, providing our perfect combination of retroflection of the endoscope tip (up to 210°) and stability of the insertion tube.

### **Enhanced visualization**

Improved visualization behind folds for increased detection of lesions during colonoscopy with an exceptional close-up view of the mucosa.

# Ergonomically designed to improve operator and patient comfort

Enhanced handling capabilities during long and complex colonoscopy procedures leading to enriched quality outcomes and patient satisfaction.

## RetroView™ Video Colonoscopes EC34-i10T

Туре	EC34-i10TL	EC34-i10TF	EC34-i10TM
Direction of view	Forward	Forward	Forward
Field of view [°]	140	140	140
Depth of field [mm]	2–100	2–100	2–100
Tip angulation up-down [°]	210–180	210–180	210–180
Tip deflection right-left [°]	160–160	160–160	160–160
Distal end width [mm]	10.5	10.5	10.5
Insertion tube width [mm]	11.6	11.6	11.6
Minimum instrument channel width [mm]	3.2	3.2	3.2
Image Sensor	HD+	HD+	HD+
Insertion tube working length [mm]	1,700	1,500	1,300

EC REP EMEA Headquarter Germany

PENTAX Europe GmbH

Julius-Vosseler-Straße 104 . 22527 Hamburg Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13

E-mail: info.emea@pentaxmedical.com . www.pentaxmedical.com

